Figure 6


Figure 7


Figure 8



Figure 9

Rayner intraocular lenses (IOLs) are single piece optical devices, manufactured from Rayacryl (hydroxyethyl methacrylate/methyl methacrylate copolymer with UV blocker). No substances are released from the material into the eye. These devices are designed to be surgically implanted in the human eye as a replacement for the crystalline lens, and are intended for placement in the capsular bag following either phacoemulsification or manual extracapsular cataract extraction (ECCE).

Rayner IOLs are intended to provide adjustment to the dioptric power of the eye. Additionally, aspheric models are aberration neutral and therefore do not add to the spherical aberration of the eye. Toric models are intended to provide adjustment to the astigmatism of the eye. Trifocal models are intended to provide pseudo accommodation to the eye.

IOL Material (Rayacryl) Characteristics

- Water content: 26% in equilibrium
- Refractive index: 1.46
- UV light transmission: shown in Figure 8 (UV 10% cut-off is 380 nm)
- Nd: YAG laser compatible

Indications

Rayner IOLs are indicated for the visual correction of aphakia in patients in whom a cataractous lens has been removed by phacoemulsification or extracapsular cataract extraction. These devices are intended to be placed in the capsular bag. Cataract surgery is performed by a suitably qualified surgeon in a sterile environment such as a hospital or Eye Clinic.

Contraindications

Apart from non-specific contraindications related to any form of ocular surgery, the following non-exhaustive list of specific contraindications must be respected:

- Microphthalmia
- Active ocular diseases (e.g. chronic severe uveitis, proliferative diabetic retinopathy, chronic glaucoma not responsive to medication)
- Corneal decompensation
- Endothelial insufficiency
- Children under the age of 18 years (toric and trifocal IOLs only)

Additionally for trifocal IOLs:

- Pseudo exfoliation
- Pars planitis
- Patients with ocular disorders, other than cataract, that could potentially cause future acuity losses to a level of 20/30 or worse in either eye
- Patients who are expected to require retinal laser treatment
- Inability to achieve secure placement in the designated location e.g. absence of secure peripheral anterior capsule, zonule laxity or dehiscence; posterior synechiae to the capsular bag
- Patients unlikely to adapt to simultaneous multiple retinal images
- Corneal oedema
- Greater than 1 Dioptre of pre-operative corneal astigmatism (non-toric trifocal model only)

Adverse Events

Cataract surgery for IOL implantation presents risks which the surgeon must evaluate. Potential complications of cataract surgery are;

- Secondary glaucoma
- IOL replacement or extraction
- IOL calcification
- Reduced vision
- Vitreous herniation
- Excessive intraoperative vitreous loss
- IOL decentration
- Secondary membrane
- Expulsive haemorrhage
- IOL dislocation and subluxation
- Retrolenticular membrane
- Corneal oedema
- Endophthalmitis and panophthalmitis
- Retinal detachment
- Corneal dystrophy
- Haemorrhage
- Iris atrophy
- Pupillary block
- Cystoid macular oedema
- Severe ametropia and aniseikonia
- Iridocyclitis and hyalitis
- Deviation from target refraction
- Fibrin reaction.

Warnings

- RayOne cannot be reused, as it is not designed to perform as intended after the first and only usage. Changes in mechanical, physical or chemical characteristics, under conditions of repeated use, cleaning and resterilisation, will compromise the integrity of the device.

- Healthcare professionals considering lens implantation under any of the following circumstances should weigh the potential risk/benefit ratio:
 - Recurrent ocular disease (e.g. uveitis, diabetic retinopathy, glaucoma, corneal decompensation)
 - Previous ocular surgery
 - Vitreous loss
 - Iris atrophy
 - Severe aniseikonia
 - Ocular haemorrhage
 - Macular degeneration
 - Zonular dehiscence (for patients at risk of zonular dehiscence, it is recommended that a capsular tension ring (CTR) is inserted to support the capsular bag).
 - Ruptured posterior capsule
 - Patients in whom the intraocular lens may affect the ability to observe, diagnose or treat posterior segment diseases.
 - Surgical difficulties at the time of cataract extraction which might increase the potential for complications (e.g. persistent bleeding, significant iris damage, uncontrolled positive pressure or significant vitreous prolapse or loss).
 - A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible.
 - Circumstances that would result in damage to the endothelium during implantation.
 - Suspected microbial infection
- Non-toothed, polished instruments must be used if handling the IOL.
- Do not allow the IOL to dehydrate during the procedure.
- Do not use balanced salt solution (saline) as the sole lubricating agent, always use in combination with an ophthalmic viscosurgical device (OVD).
- Do not attempt to disassemble, modify or alter this device or any of its components, as this can significantly affect the function and/or structural integrity of the design.
- In the rare event that lens damage is observed post injection, the clinician must make an assessment of the severity as relates to likely optical outcomes and take corrective action if appropriate (i.e., replacing the IOL intraoperatively may be appropriate if the IOL is significantly decentred within the capsular bag).

Precautions

- The sterility of the contents is guaranteed only if the outer tray has not been opened or damaged.
- Do not use if the pack has been damaged.
- Do not store in direct sunlight.
- Do not store the pack outside of the recommended storage conditions; store between 5°C to 25°C.
- Do not store below 20% Relative Humidity.
- Do not use after the expiry date.
- Do not attempt to reuse this device.
- Do not resterilise.
- Do not use system at temperature below 21°C (equilibration takes approximately 90 minutes from an initial temperature of 0°C).

Sterilisation and Packaging

The IOL is supplied sterile and preloaded in the RayOne delivery system within a sterilised blister pack. The sterilised blister pack is steam sterilised and should only be opened under sterile conditions.

Patient Implant Card

An implant card is included in the pack to record all implant information (the supplied labels may be used). It shall be given to the patient, with the instruction to keep this card. The card should be shown to any eye care professional the patient visits in future.

Directions for Use

The RayOne preloaded injector should only be used for the placement of Rayner IOLs into the eye.

To ensure optimum injection performance of the IOLs, the blister pack should be allowed to equilibrate to a temperature of 21°C or above before use. The use of a sodium hyaluronate based OVD is recommended. Insert the IOL into the eye within 3 minutes of closing the cartridge.

Use of RayOne

Fig 1 Completely peel back the lid of the first blister pack.

Fig 2 Carefully peel back the lid half way down the second blister pack.

Fig 3 Carefully drain the saline from the blister tray and peel off remaining lid.

Do not remove the injector from the blister tray. Do not wait more than 3 minutes before adding OVD - Dehydration risk.

Fig 4 Insert the viscoelastic cannula into the opening marked with an arrow on the cartridge and apply sufficient OVD to completely fill the cartridge.

NOTE: Rayner's patented "Lock & Roll" injector technology uses the mechanical force of OVD inserted into the indicated port to move the IOL from 'storage' to 'roll ready' position at the bottom of the cartridge; do not add OVD directly via the nozzle tip.

Fig 5 Keep the injector in the tray and close the cartridge firmly together by pushing the moving half of the cartridge (Labelled 2) against the fixed half until you hear it click closed. Check both clips have "clicked" shut and secured the cartridge.

Fig 6 Gently lift out the injector from the tray.

Fig 7 Press the plunger in a slow and controlled manner. If excessive resistance is felt this could indicate a blockage; discontinue use of the product and return the product and all packaging to Rayner. In the case of IOL rotation during ejection from the nozzle, gently rotate the injector in the opposite direction to counteract any movement. Stop depressing the plunger when the IOL exits the nozzle. Discard the injector after use. Following implantation, irrigate/aspirate to eliminate any OVD residues from the eye, especially behind the IOL.

IOL Placement

Ensure that the IOL is implanted in the capsular bag with the correct anterior/posterior placement. Correct anterior/posterior placement is achieved when the haptics sweep away from the optic in a counter-clockwise direction (anterior view Figure 10). The IOL can be dialed clockwise as is usual for posterior chamber IOLs. The anterior continuous curvilinear capsulorhexis should be positioned just covering the 360° anterior edge of the IOL optic by 0.5 to 1.0 mm. Correct rotational orientation for toric IOLs is achieved when the IOL's axis marks (the lowest IOL power meridian) aligns with the cornea's highest power meridian. The IOL axis marks are therefore oriented in parallel to (or on) the steeper corneal meridian.

Calculation of IOL Power

The surgeon should preoperatively determine the power of the IOL to be implanted. Lens power calculation methods are described in the following references:

1. Retzlaff J., Sanders D. & Kraff M. Lens Implant Power Calculation, A Manual for Ophthalmologists & Biometrists - Third Edition, 1990.
2. Holladay J. A Three-part System For Refining Intraocular Lens Power Calculations. J. Cataract Refract. Surg. V14:17-24, 1988.
3. Holladay J. Standardizing Constants For Ultrasonic Biometry, Keratometry & IOL Power Calculations. J. Cataract Refract. Surg. V23:1356-1370, 1997.
4. Hoffer K. The Hoffer Q Formula: A Comparison of Theoretic and Regression Formulas. J. Cataract Refract. Surg. V19:700-712, 1993; ERRATA, 20:677, 1994; ERRATA, 33:2, 2007

Method of Disposal

Products contaminated with biological fluids shall be disposed of in accordance with local guidance relating to medical waste. Packaging should be disposed of in accordance with local guidance relating to recycling.

Electromagnetic Compatibility

Rayner IOLs are considered safe for magnetic resonance imaging (MRI), and do not pose an increased risk during diagnostic investigation or therapeutic treatment.

Figure 10

IOL Haptic Orientation

Symbol/Explanation

	Single use, do not reuse		Do not re-sterilise
	Steam sterilised		Do not use if the sterile barrier system or packaging is damaged
	Store between 5-25°C temperature		Keep away from sunlight
	Keep dry		Consult Instructions for use
	Medical Device		European Authorized Representative

Summary of Safety and Clinical Performance

Rayner has prepared a Summary of Safety and Clinical Performance (SSCP) document which has been validated in English by BSI Netherlands (Notified Body number 2797). BSI will make this document available on the Eudamed database which can be found at ec.europa.eu/tools/Eudamed. The SSCP can be found using the Basic UDI-DI references:

Product Name	Model Number	Basic UDI-DI
RayOne Aspheric	RA0600C	05029867004070
RayOne EMV	RA0200E	05029867004254
RayOne Spheric	RA0100C	05029867004087
RayOne Toric	RA0610T	05029867004094
RayOne Trifocal	RA0603F	05029867004100
RayOne Trifocal Toric	RA0613Z	05029867004117
RayOne EMV Toric	RA0210T	05029867003912

Customer Feedback

Any serious incidents shall be immediately reported to Rayner and the competent authority of the Member State in which the user is established. For competent authority information, visit www.rayner.com/incident.

Please report any feedback, adverse events/complaints to:
Rayner, 10 Dominion Way, Worthing, West Sussex, BN14 8AQ, United Kingdom.

Tel: +44 (0) 1903 258900

Fax: +44 (0) 1903 258901

Email: feedback@rayner.com

A product lifetime of twenty years has been assigned to these devices, corresponding to the shelf-life of the product, typical age of the patient and typical life expectancy.

The intraocular lens in this pack is for implantation into the capsular bag. A non-exhaustive list of risks, instructions for use, cautions and precautions and warnings are provided in this IFU which should be given to the patient prior to surgery. Rayner hereby disclaims all liability associated with any injuries or harm that may be suffered by the patient caused by: The method or technique of implantation used by the surgeon if the manufacturer's recommendations and instructions were not followed; an unsuitable prescription, selection or implantation of an unsuitable intraocular lens.

©2023 Rayner. Rayner, RayOne and Rayacryl are proprietary marks of Rayner. Rayner Intraocular Lenses Limited, 10 Dominion Way, Worthing, West Sussex, BN14 8AQ, United Kingdom. Registered in England: 615539.

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1. Retzlaff J., Sanders D. & Kraff M. Lens Implant Power Calculation, A Manual for Ophthalmologists & Biometrists - Third Edition, 1990. 2. Holladay J. A Three-part System For Refining Intraocular Lens Power Calculations. J. Cataract Refract. Surg. V14:17-24, 1988. 3. Holladay J. Standardizing Constants For Ultrasonic Biometry, Keratometry & IOL Power Calculations. J. Cataract Refract. Surg. V23:1356-1370, 1997. 4. Hoffer K. The Hoffer Q Formula: A Comparison of Theoretic and Regression Formulas. J. Cataract Refract. Surg. V19:700-712, 1993; ERRATA, 20:677, 1994; ERRATA, 33:2, 2007

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1. Retzlaff J., Sanders D. & Kraff M. Lens Implant Power Calculation, A Manual for Ophthalmologists & Biometrists - Third Edition, 1990
2. Holladay J. A Three-part System For Refining Intraocular Lens Power Calculations. J. Cataract Refract. Surg. V14:17-24, 1988
3. Holladay J. Standardizing Constants For Ultrasonic Biometry, Keratometry & IOL Power Calculations. J. Cataract Refract. Surg. V23:1356-1370, 1997
4. Hoffer K. The Hoffer Q Formula: A Comparison of Theoretic and Regression Formulas. J. Cataract Refract. Surg. V19:700-712, 1993; ERRATA, 20:677, 1994; ERRATA, 33:2, 2007

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