
# Ray One

HYDROPHOBIC

Preloaded Hydrophobic Acrylic  
IOL Injection System

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

# Use of Preloaded Hydrophobic Acrylic IOL Injection System


**Figure 1**


**Figure 2**


**Figure 3**


**Figure 4**


**Figure 5**

**Figure 6**


**Figure 7**

**Figure 8**

**IOL Material UV Light Transmission**

— RAO800C

■ RAO850B

Rayner RayOne Hydrophobic intraocular lenses (IOLs) are single piece optical devices, manufactured from a photopolymer composition of hydrophobic polyurethane acrylates and methacrylates with UV blocker. RayOne Hydrophobic BLF (Model RA0850B) additionally contains a blue light filter. These devices are designed to be surgically implanted in the human eye as a replacement for the crystalline lens, and are intended for placement in the capsular bag following either phacoemulsification or manual extracapsular cataract extraction (ECCE). Rayner IOLs are intended to provide adjustment to the dioptric power of the eye. Additionally, aspheric models are aberration neutral and therefore do not add to the spherical aberration of the eye.

### **IOL Material Characteristics**

- Refractive index: 1.51 (800C, 850B)
- UV light transmission: shown in Figure 8 (UV 10% cut-off is 385 nm [800C], 386 nm [850B])
- Nd:YAG laser compatible

### **Indications**

Rayner IOLs are indicated for the visual correction of aphakia in patients in whom a cataractous lens has been removed by phacoemulsification or extracapsular cataract extraction. These devices are intended to be placed in the capsular bag.

### **Contraindications**

Apart from non-specific contraindications related to any form of ocular surgery, the following non-exhaustive list of specific contraindications must be respected:

- Microphthalmia
- Active ocular diseases (e.g. chronic severe uveitis, proliferative diabetic retinopathy, chronic glaucoma not responsive to medication)
- Corneal decompensation
- Endothelial insufficiency
- Children under the age of 18

### **Adverse Events**

Cataract surgery for IOL implantation presents risks which the surgeon must evaluate. Potential complications of cataract surgery are:

- Secondary glaucoma
- IOL replacement or extraction
- Precipitates
- Reduced vision
- Vitreous herniation
- Excessive intraoperative vitreous loss
- IOL decentration
- Secondary membrane
- Expulsive haemorrhage
- IOL dislocation and subluxation
- Retrolenticular membrane
- Corneal oedema
- Endophthalmitis and panophthalmitis
- Retinal detachment
- Corneal dystrophy
- Haemorrhage
- Iris atrophy
- Pupillary block
- Cystoid macular oedema
- Severe ametropia and aniseikonia
- Iridocyclitis and hyalitis
- Deviation from target refraction
- Fibrin reaction

### **⚠ Warnings**

- RayOne Hydrophobic cannot be reused, as it is not designed to perform as intended after the first and only usage. Changes in mechanical, physical or chemical characteristics, under conditions of repeated use, cleaning and resterilisation, will compromise the integrity of the device
- Healthcare professionals considering lens implantation under any of the following circumstances should weigh the potential risk/benefit ratio:

- Recurrent ocular disease (e.g. uveitis, diabetic retinopathy, glaucoma, corneal decompensation)
- Previous ocular surgery
- Vitreous loss
- Iris atrophy
- Severe aniseikonia
- Ocular haemorrhage
- Macular degeneration
- Zonular dehiscence (for patients at risk of zonular dehiscence, it is recommended that a capsular tension ring (CTR) is inserted to support the capsular bag).
- Ruptured posterior capsule
- Patients in whom the intraocular lens may affect the ability to observe, diagnose or treat posterior segment diseases
- Surgical difficulties at the time of cataract extraction which might increase the potential for complications (e.g. persistent bleeding, significant iris damage, uncontrolled positive pressure or significant vitreous prolapse or loss)
- A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible
- Circumstances that would result in damage to the endothelium during implantation
- Suspected microbial infection

- Non-toothed, polished instruments must be used if handling the IOL
- Do not use unsterile surgical instruments or instruments that may carry a risk of contamination
- Do not use balanced salt solution (saline) as the sole lubricating agent, always use in combination with an ophthalmic viscosurgical device (OVD)
- Do not attempt to disassemble, modify or alter this device or any of its components, as this can significantly affect the function and/or structural integrity of the design
- Do not attempt to re-open and re-close the flaps after closing the cartridge for the first time.

## Precautions

- The sterility of the contents is guaranteed only if the outer tray has not been opened or damaged
- Do not use if the pack has been damaged
- Do not store in direct sunlight
- Do not store the pack outside of the recommended storage conditions; store between 5°C to 25°C
- Do not use after the expiry date
- Do not attempt to reuse this device
- Do not re-sterilise
- Do not use system at temperature below 18°C (allow system to equilibrate to room temperature)

## Sterilisation and Packaging

The IOL is supplied sterile and preloaded in the RayOne delivery system within a sterilised blister pack. The blister pack is steam sterilised and should only be opened under sterile conditions. An implant card is included in the pack to record all implant information (the supplied labels may be used). It shall be given to the patient, with the instruction to keep this card. The card should be shown to any eye care professional the patient visits in future.

## Directions for Use

The RayOne preloaded injector should only be used for the placement of Rayner IOLs into the eye.

To ensure optimum injection performance of the IOLs, the blister pack should be allowed to equilibrate to a temperature of 18°C or above before use. The use of a sodium hyaluronate based OVD is recommended.

## Use of RayOne

**Fig 1** Completely peel back the lid of the first blister pack.

**Fig 2** Carefully peel back the lid of the second blister pack.

**Fig 3** Carefully drain the saline from the blister pack and peel off remaining lid. Do not remove the injector from the blister tray.

**Fig 4** Do not remove the injector from the blister tray. Insert the viscoelastic cannula into the opening marked with an arrow on the cartridge and apply sufficient OVD to completely fill the cartridge.

NOTE: Rayner's patented "Lock & Roll" injector technology uses the mechanical force of OVD inserted into the indicated port to move the IOL from 'storage' to 'roll ready' position at the bottom of the cartridge; do not add OVD directly via the nozzle tip.

**Fig 5** Keep the injector in the tray and close the cartridge firmly together by pushing the moving half of the cartridge (Labelled 2) against the fixed half until you hear it click closed. Check both clips have "clicked" shut and secured the cartridge.

**Fig 6** Gently lift out the injector from the tray.

**Fig 7** Press the plunger in a slow and controlled manner. If excessive resistance is felt this could indicate a blockage; discontinue use of the product and return the product and all packaging to Rayner. In the case of IOL rotation during ejection from the nozzle, gently rotate the injector in the opposite direction to counteract any movement. Stop depressing the plunger when the IOL exits the nozzle. Discard the injector after use.

Following implantation, irrigate/aspirate to eliminate any OVD residues from the eye, especially behind the IOL.

## IOL Placement

Ensure that the IOL is implanted in the capsular bag with the correct anterior/posterior placement. Correct anterior/posterior placement is achieved when the haptics sweep away from the optic in a counter-clockwise direction (anterior view Figure 9). The IOL can be dialled clockwise as is usual for posterior chamber IOLs. The anterior continuous curvilinear capsulorhexis should be positioned just covering the 360° anterior edge of the IOL optic by 0.5 to 1.0 mm.

**Figure 9**

### IOL Haptic Orientation

## Calculation of IOL Power

The surgeon should preoperatively determine the power of the IOL to be implanted. Lens power calculation methods are described in the following references:

1. Retzlaff J., Sanders D. & Kraff M. Lens Implant Power Calculation, A Manual for Ophthalmologists & Biometrists - Third Edition, 1990.
2. Holladay J. A Three-part System For Refining Intraocular Lens Power Calculations. J. Cataract Refract. Surg. V14:17-24, 1988.
3. Holladay J. Standardizing Constants For Ultrasonic Biometry, Keratometry & IOL Power Calculations. J. Cataract Refract. Surg. V23:1356-1370, 1997.
4. Hoffer K. The Hoffer Q Formula: A Comparison of Theoretic and Regression Formulas. J. Cataract Refract. Surg. V19:700-712, 1993; ERRATA, 20:677, 1994; ERRATA, 33:2, 2007

## Method of Disposal

Products contaminated with biological fluids shall be disposed of in accordance with local guidance relating to medical waste. Packaging should be disposed of in accordance with local guidance relating to recycling.

## Electromagnetic Compatibility

Rayner IOLs are considered safe for magnetic resonance imaging (MRI), and do not pose an increased risk during diagnostic investigation or therapeutic treatment.

## Symbol/Explanation

	Single use, do not reuse		Do not re-sterilise
	Steam Sterilised		Do not use if the sterile barrier system or packaging is damaged
	Store between 5-25°C temperature		Keep away from sunlight
	Keep dry		Consult Instructions for use
	Batch code		Serial Number
	Manufacturer		Use-by date
	Medical Device		European Authorized Representative

## Customer Feedback

Any serious incidents shall be immediately reported to Rayner and the competent authority of the Member State in which the user is established.

For competent authority information, visit [www.rayner.com/incident](http://www.rayner.com/incident).

Please report any feedback, adverse events/complaints to:

Rayner, 10 Dominion Way, Worthing, West Sussex, BN14 8AQ, United Kingdom.

**Tel: +44 (0) 1903 258900**

**Fax: +44 (0) 1903 258901**

**Email: [feedback@rayner.com](mailto:feedback@rayner.com)**

A product lifetime of twenty years has been assigned to these devices, corresponding to the shelf-life of the product, typical age of the patient and typical life expectancy.

The intraocular lens in this pack is for implantation into the capsular bag.

A non-exhaustive list of risks, instructions for use, cautions and precautions and warnings are provided in this IFU which should be given to the patient prior to surgery. Rayner hereby disclaims all liability associated with any injuries or harm that may be suffered by the patient caused by: The method or technique of implantation used by the surgeon if the manufacturer's recommendations and instructions were not followed; an unsuitable prescription, selection or implantation of an unsuitable intraocular lens.

©2021 Rayner. Rayner and RayOne are proprietary marks of Rayner.

Rayner Intraocular Lenses Limited, 10 Dominion Way, Worthing,

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
1. Retzlaff J., Sanders D. & Kraff M. Lens Implant Power Calculation, A Manual for Ophthalmologists & Biometrists - Third Edition, 1990. 2. Holladay J. A Three-part System For Refining Intraocular Lens Power Calculations. J. Cataract Refract. Surg. V14:17-24, 1988. 3. Holladay J. Standardizing Constants For Ultrasonic Biometry, Keratometry & IOL Power Calculations. J. Cataract Refract. Surg. V23:1356-1370, 1997. 4. Hoffer K. The Hoffer Q Formula: A Comparison of Theoretic and Regression Formulas. J. Cataract Refract. Surg. V19:700-712, 1993; ERRATA, 20:677, 1994; ERRATA, 33:2, 2007

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1. Retzlaff J., Sanders D. & Kraff M. Lens Implant Power Calculation, A Manual for Ophthalmologists & Biometrists - Third Edition, 1990
2. Holladay J. A Three-part System For Refining Intraocular Lens Power Calculations. J. Cataract Refract. Surg. V14:17-24, 1988
3. Holladay J. Standardizing Constants For Ultrasonic Biometry, Keratometry & IOL Power Calculations. J. Cataract Refract. Surg. V23:1356-1370, 1997
4. Hoffer K. The Hoffer Q Formula: A Comparison of Theoretic and Regression Formulas. J. Cataract Refract. Surg. V19:700-712, 1993; ERRATA, 20:677, 1994; ERRATA, 33:2, 2007

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
1. Retzlaff J., Sanders D. & Kraff M. Lens Implant Power Calculation, A Manual for Ophthalmologists & Biometrists – Third Edition, 1990
2. Holladay J. A Three-part System For Refining Intraocular Lens Power Calculations. J. Cataract Refract. Surg. V14:17-24, 1988
3. Holladay J. Standardizing Constants For Ultrasonic Biometry, Keratometry & IOL Power Calculations. J. Cataract Refract. Surg. V23:1356-1370, 1997
4. Hoffer K. The Hoffer Q Formula: A Comparison of Theoretic and Regression Formulas. J. Cataract Refract. Surg. V19:700-712, 1993; ERRATA, 20:677, 1994; ERRATA, 33:2, 2007

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