# INTRAOCULAR LENS MODELS

*C-flex*

*Superflex*


*C-flex Aspheric*

*Superflex Aspheric*

*T-flex Aspheric*

*M-flex*

*M-flex T*






Intraocular lenses  
with single use injectors

CE  
2797

 Rayner

# Use of injector




Figure 1. M-flex +4 Dioptre add through focus MTF



Figure 2. U V light transmission

## Description

This Pack contains:

- One injectable intraocular lens
- One single use, soft-tipped injector

The injector provided in this pack shall only be used in combination with the lens provided in this pack.

Rayner intraocular lenses (IOLs) are single piece optical devices, manufactured from Rayacryl (hydroxyethyl methacrylate/methyl methacrylate copolymer with UV blocker). These devices are designed to be surgically implanted in the human eye as a replacement for the crystalline lens, and are intended for placement in the capsular bag following either phacoemulsification or manual extracapsular cataract extraction (ECCE).

Rayner IOLs are intended to provide adjustment to the dioptric power of the eye. Additionally, aspheric models are aberration neutral and therefore do not add to the spherical aberration of the eye; toric models are intended to provide adjustment to astigmatism of the eye; multifocal models are intended to provide pseudo accommodation to the eye.

**Note:** For through focus MTF performance of multifocal models, refer to Figure 1. Through focus MTF performance of a +4 Dioptre add *M-flex* IOL was measured in a standardised model eye (ISO 11979) at 50 lp/mm as the focus shifted from far objects (0 Dioptres) to increasingly nearer objects (higher MTF indicates better optical performance).

## IOL Material (Rayacryl) Characteristics

- Water content: 26% in equilibrium
- Refractive index: 1.46
- UV light transmission: shown in Figure 2 (UV 10% cut-off is 380nm)
- Nd:YAG laser compatible

## Indications

- Aphakia following cataract surgery
- Senile cataract
- Traumatic cataract
- Congenital or juvenile cataract (spherical and aspheric models only)
- Correction of regular corneal astigmatism >0.75 Dioptres (toric models only)

## Contraindications

Apart from non-specific contraindications related to any form of ocular surgery, the following non-exhaustive list of specific contraindications must be respected:

- Microphthalmia
- Active ocular diseases (chronic severe uveitis, proliferative diabetic retinopathy, chronic glaucoma not responsive to medication)
- Corneal decompensation or endothelial insufficiency
- Children under the age of 18 years (toric and multifocal IOLs only)

### **Additionally for multifocal IOLs:**

- Pseudo exfoliation
- High myopia
- Pars planitis
- Patients with ocular disorders, other than cataract, that could potentially cause future acuity losses to a level of 20/30 or worse in either eye
- Patients who are expected to require retinal laser treatment
- Inability to achieve secure placement in the designated location e.g. absence of a secure peripheral anterior capsule, zonule laxity or dehiscence; posterior synechiae to the capsular bag
- Patients unlikely to adapt to simultaneous multiple retinal images
- Greater than 1 Dioptre of pre-operative corneal astigmatism (non-toric multifocal model only)

## Complications Related to IOL Surgery

A non-exhaustive list of intraoperative & postoperative complications

- Secondary glaucoma
- IOL replacement or extraction
- IOL precipitates
- Reduced vision
- Vitreous herniation
- Excessive intraoperative vitreous loss
- IOL decentration
- Secondary membrane
- Expulsive haemorrhage
- IOL dislocation and subluxation
- Retrolenticular membrane
- Corneal oedema
- Endophthalmitis and panophthalmitis
- Retinal detachment
- Corneal dystrophy
- Haemorrhage
- Iris atrophy
- Pupillary block
- Cystoid macular oedema
- Severe ametropia and aniseikonia
- Iridocyclitis and hyalitis
- Deviation from target refraction
- Fibrin reaction

## Warnings

- Single use IOLs and single use injectors cannot be reused, as they are not designed to perform as intended after the first and only usage. Changes in mechanical, physical or chemical characteristics, under conditions of repeated use, cleaning and resterilisation, will compromise the integrity of IOLs and injectors.
- Do not use after the expiry date.
- Do not resterilise.
- Non-toothed, polished instruments must be used if handling the IOL.
- Do not allow the IOL to contact substances that are unsterile or ocular-incompatible prior to placement into the eye. Do not use unsterile surgical instruments or instruments that may carry a risk of contamination.
- Do not allow the IOL to dehydrate during the procedure.
- Do not use balanced salt solution (saline) as the sole lubricating agent during injection, but always in combination with an ophthalmic viscosurgical device (OVD).
- Once closed, do not reopen the flaps of the injector.
- Following implantation, irrigate/aspirate to eliminate any OVD residues from the eye, especially behind the IOL.

## Precautions

A risk/benefit analysis must be performed before confirming a patient as a candidate for a Rayner IOL, if they are suffering from any of the following conditions:

- Recurrent ocular disease (e.g. uveitis, diabetic retinopathy, glaucoma, corneal decompensation)
- Previous ocular surgery
- Non-age related cataract
- Vitreous loss
- Iris atrophy
- Severe Aniseikonia
- Ocular Haemorrhage
- Macular degeneration
- Zonular dehiscence (for patients at risk of zonular dehiscence, it is recommended that a capsular tension ring (CTR) is inserted to support the capsular bag)
- Ruptured posterior capsule
- Patients in whom the intraocular lens may affect the ability to observe, diagnose, or treat posterior segment diseases
- Surgical difficulties at the time of cataract extraction which might increase the potential for complications (e.g. persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss)
- A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible
- Circumstances that would result in damage to the endothelium during implantation
- Suspected microbial infection

## Sterilisation and Packaging

The IOL is supplied sterile in a sealed blister pack containing 0.9% saline solution. The sterilised blister pack is steam sterilised and should only be opened under sterile conditions. An implant card is included in the pack to record all implant information (the supplied labels may be used). It shall be given to the patient, with the instruction to keep this card. The card should be shown to any eye care professional the patient visits in future.

## Directions for Use

Rayner injectors should only be used for the placement of Rayacryl IOLs into the eye. To ensure optimum injection performance of the IOLs, the blister pack should be allowed to equilibrate to a temperature of 21°C or above before use (equilibration takes approximately 30 minutes from an initial temperature of 0°C). Load the injector under well lit conditions or using magnification. The use of a sodium hyaluronate-based OVD is recommended. Load the lens into the injector immediately after removal from the blister pack and insert into the eye within 3 minutes of loading. For loading details, refer to the 'Use of Injector' section of this IFU.

## IOL Placement

The surgeon should ensure that the IOL is implanted in the capsular bag, with the correct anterior/posterior placement. Correct anterior/posterior placement is achieved when the haptics sweep away from the optic in a counter-clockwise direction (anterior view Figure 10). This means that the IOL can most easily be dialled clockwise as is usual for posterior chamber IOLs. The anterior continuous curvilinear capsulorhexis should be positioned just covering the 360° anterior edge of the IOL optic by 0.5 to 1.0 mm. Correct rotational orientation for toric IOLs is achieved when the IOL's axis marks (the lowest IOL power meridian) aligns with the cornea's highest power meridian. The IOL axis marks are therefore oriented in parallel to (or on) the steeper corneal meridian.

**Note:** The toric axis marks are on the anterior IOL surface.

## Use of Injector (Refer to Diagrams 1 - 8)

**Step 1** Aseptically transfer the injector to the sterile field by tipping it from the tray. Fully retract the plunger ensuring that the soft tip does not protrude into the loading bay. Do not advance the plunger prior to loading the lens.

**Step 2** Open the loading bay flap fully to 90° and apply OVD inside the nozzle and to both grooves of the loading bay.

**Step 3** Carefully peel back the foil lid of the lens blister. Gently lift out the lens using parallel tipped, non-serrated forceps e.g. Kelman-McPherson. Rinse the lens with sterile balanced salt solution.

**Step 4** Position the lens centrally in the loading bay in a **"reverse - S" configuration**. Ensure that the nearest edge of the optic is securely under the edge (lip) of the loading bay as shown.

**Step 5** Hold open the flap and gently press down on the lens with closed forceps to ensure the furthest edge of the optic is securely under the edge (lip) of the flap as shown (A). Ensure the haptics are also securely under the edges within the loading bay (B and C) and orientated toward the optic.

**Step 6** **Ensure that NO parts of the optic or haptics of the lens are projecting outside the edges of the loading bay.** While keeping the lens in position with open forceps and placing gentle down pressure on the optic, carefully close the flaps of the injector locking them firmly together. Any resistance could indicate a trapped lens.

**Step 7** Ensure that no parts of the optic or haptics are trapped between the flaps. Advance the plunger in a slow controlled manner. Anticipate an initial slight resistance. Excessive resistance could indicate a trapped lens. If excessive resistance is felt, fully retract the plunger and then advance until in contact with the lens again. If the lens causes a blockage in the injector system, discard the injector.

**Step 8** Continue the injection in a **slow and controlled manner**. Do not exert excessive force on the plunger. When the lens exits the nozzle, stop depressing the plunger and do not retract the plunger. Discard the injector after use.

**Figure 10**

### IOL Haptic Orientation


## Calculation of IOL Power

The surgeon should preoperatively determine the power of the IOL to be implanted. This can be calculated variously from the corneal radius, the depth of the anterior chamber and the axial length of the eye according to formulae described in the following references:

1. Retzlaff J., Sanders D. & Kraff M. Lens Implant Power Calculation, A Manual for Ophthalmologists & Biometrists - Third Edition, 1990
2. Holladay J. A Three-part System For Refining Intraocular Lens Power Calculations. J. Cataract Refract. Surg. V14:17-24, 1988
3. Holladay J. Standardizing Constants For Ultrasonic Biometry, Keratometry & IOL Power Calculations. J. Cataract Refract. Surg. V23:1356-1370, 1997
4. Hoffer K. The Hoffer Q Formula: A Comparison of Theoretic and Regression Formulas. J. Cataract Refract. Surg. V19:700-712, 1993; ERRATA, 20:677, 1994; ERRATA, 33:2, 2007

## Symbol/Explanation

	Consult Instructions for Use		Caution Consult the Instructions for Use for warnings and precautions
	Single use, do not reuse		Do not re-sterilise
	Steam sterilised (IOL)		Ethylene Oxide sterilised (injector)
	Store between 0-45°C		Do not use if the sterile barrier system or packaging is damaged
	Keep dry		Keep away from sunlight

**CE**  
**2797**

## Customer Feedback

Please report any feedback, adverse events/complaints to:  
Rayner Intraocular Lenses Limited, The Ridley Innovation Centre,  
10 Dominion Way, Worthing, West Sussex, BN14 8AQ, United Kingdom.

**Tel: +44 (0) 1903 258900**

**Fax: +44 (0) 1903 258901**

**Email: [feedback@rayner.com](mailto:feedback@rayner.com)**

The intraocular lens in this pack is for implantation into the capsular bag. A non-exhaustive list of risks, instructions for use, cautions and precautions and warnings are provided in this IFU which should be given to the patient prior to surgery. Rayner hereby disclaims all liability associated with any injuries or harm that may be suffered by the patient caused by: The method or technique of implantation used by the surgeon if the manufacturer's recommendations and instructions were not followed; an unsuitable prescription, selection or implantation of an unsuitable intraocular lens.

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
1. Retzlaff J., Sanders D. & Kraff M. Lens Implant Power Calculation, A Manual for Ophthalmologists & Biometrists - Third Edition,1990
2. Holladay J. A Three-part System For Refining Intraocular Lens Power Calculations. J. Cataract Refract. Surg. V14:17-24,1988
3. Holladay J. Standardizing Constants For Ultrasonic Biometry, Keratometry & IOL Power Calculations. J. Cataract Refract. Surg. V23:1356-1370,1997
4. Hoffer K. The Hoffer Q Formula: A Comparison of Theoretic and Regression Formulas. J. Cataract Refract. Surg. V19:700-712, 1993; ERRATA, 20:677, 1994; ERRATA, 33:2,2007

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**CE**  
**2797**

## Feedback del cliente

Comunicare feedback, eventuali problemi / lamentele a:  
Rayner Intraocular Lenses Limited, The Ridley Innovation Centre,  
10 Dominion Way, Worthing, West Sussex, BN14 8AQ, United Kingdom.

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1. Retzlaff J., Sanders D. & Kraff M. Lens Implant Power Calculation, A Manual for Ophthalmologists & Biometrists - Third Edition, 1990
2. Holladay J. A Three-part System For Refining Intraocular Lens Power Calculations. J. Cataract Refract. Surg. V14:17–24, 1988
3. Holladay J. Standardizing Constants For Ultrasonic Biometry, Keratometry & IOL Power Calculations. J. Cataract Refract. Surg. V23:1356–1370, 1997
4. Hoffer K. The Hoffer Q Formula: A Comparison of Theoretic and Regression Formulas. J. Cataract Refract. Surg. V19:700–712, 1993; ERRATA, 20:677, 1994; ERRATA, 33:2, 2007

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Rayner Intraocular Lenses Limited, The Ridley Innovation Centre,  
10 Dominion Way, Worthing, West Sussex, BN14 8AQ, Wielka Brytania.

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1. Retzlaff J., Sanders D. & Kraff M. Lens Implant Power Calculation, A manual for Ophthalmologists & Biometrists - Third Edition, 1990
2. Holladay J. A Three-part System For Refining Intraocular Lens Power Calculations. J. Cataract Refract. Surg. V14:17-24, 1998
3. Holladay J. Standardizing Constants For Ultrasonic Biometry, Keratometry & IOL Power Calculations. J. Cataract Refract. Surg. V23:1356-1370, 1997
4. Hoffer K. The Hoffer Q Formula: A Comparison of Theoretic and Regression Formulas. J. Cataract Refract. Surg. V19:700-712, 1993.

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©2019 Rayner. Rayner и Rayacryl са фирмени марки на фирма Rayner. Rayner Intraocular Lenses Limited, 10 Dominion Way, Worthing, West Sussex, BN14 8AQ, United Kingdom. Регистрирана в Англия: 615539.