


# RayOne<sup>®</sup> Spheric Technical Information



# RayOne® Spheric Technical Information

Model Name	RayOne® Spheric
Model Number	RAO100C
Power Range	-10.0 D to +7.0 D (1.0 D increments, inc. plano) +8.0 D to +30.0 D (0.5 D increments) +31.0 D to +34.0 D (1.0 D increments)
Delivery System Type	Fully preloaded IOL injection system
Incision Size	Sub 2.2 mm

Delivery System	
Injector Type	Single use, fully preloaded IOL injection system
Nozzle Size	1.65 mm
Bevel Angle	45°
Lens Delivery	Single handed plunger

Spheric Monofocal IOL	
Material	Single piece Rayacryl® hydrophilic acrylic
Water Content	26% in equilibrium
UV Protection	Benzophenone UV absorbing agent
UV Light Transmission	UV 10% cut-off is 380 nm
Refractive Index	1.46
ABBE	56
Overall Diameter	12.50 mm
Optic Diameter	6.00 mm
Optic Shape	Biconvex (positive powers), Biconcave (negative powers)
Optic Edge Design	Amon-Apple 360° enhanced square edge
Haptic Angulation	0°, uniplanar
Haptic style	Closed loop with anti-vaulting haptic (AVH) technology

Estimated Constants for Optical Biometry					
SRK/T	Haigis			HofferQ	Holladay
A-constant	a0	a1	a2	pACD	SF
118.6	1.17	0.40	0.10	5.32	1.56

For Contact Ultrasound, the estimated A-constant is 118.0

Please note that the constants indicated for all Rayner lenses are estimates and are for guidance purposes only. Surgeons must always expect to personalise their own constants based on initial patient outcomes, with further personalisation as the number of eyes increases.

Discover why RayOne® is in a class of its own visit [rayner.com/rayone](http://rayner.com/rayone)



Regulatory restrictions may apply in certain markets.  
©Rayner 2018. RayOne is a registered mark of Rayner.  
EC 2018-133 12/18

GLOBAL HQ & UK OFFICE  
RAYNER INTRAOCULAR LENSES LIMITED  
10 Dominion Way  
Worthing, West Sussex  
BN14 8AQ, United Kingdom

Tel: +44 (0) 1903 258 900  
Email: [iol\\_enquiries@rayner.com](mailto:iol_enquiries@rayner.com)  
Web: [rayner.com](http://rayner.com)

UK CUSTOMER SERVICE  
Tel: +44 (0) 1903 854 560  
Email: [orders@rayner.com](mailto:orders@rayner.com)

GERMANY CUSTOMER SERVICE  
Tel: +49 30 6290783 - 0  
Email: [info.de@rayner.com](mailto:info.de@rayner.com)

ITALY CUSTOMER SERVICE  
Tel: +39 800 972 485  
Email: [ordini@rayner.com](mailto:ordini@rayner.com)

SPAIN CUSTOMER SERVICE  
Tel: +349 1048 0055  
Email: [pedidos@rayner.com](mailto:pedidos@rayner.com)

PORTUGAL CUSTOMER SERVICE  
Tel: +351 221 206 400  
Email: [pedidos@rayner.com](mailto:pedidos@rayner.com)