

Optimise cost and quality of

Cataract Surgery

Premium Single-Use Recyclable Instruments



Grateful for the support from surgeons worldwide



Dr. Vincent Qin, MD, MBA, MPH, FEBO Ophthalmology & Eye Surgery Consultant

"Since using the HASA OPTIX Cataract Sets, I never have to worry about missing or broken instruments. I can fully perform cataract surgery in the best material conditions."



Dr. Eric DonnenfeldClinical Professor of Ophthalmology NYU at OCLI


"These instruments are basically an extension of my hands, in the patient's eye, and I feel a level of comfort that I haven't felt otherwise."



Dr. Ivo Guber, MD

Owner & Medical Director bei augenchirurgie.ch Winterthur, Zurich, Switzerland

"Excellent quality of instruments for unbeatable price."



Ms. Masara Laginaf

Consultant & Ophthalmic Surgeon Specialized in Cornea, cataract and Refractive Surgery at OCL Vision London (UK)

"Created with a strong ethical position to improve sustainabilty in the field of ophthalmology whilst providing streamlined solutions and beautifully engineered instruments."

HCPs providing testimonials may be speakers or paid consultants of Rayner.

A genuine desire to empower all those involved in ophthalmic surgery

Hospitals, surgical centres, and ophthalmic clinics face the challenge of working profitably whilst maintaining high levels of patient care.

HASA OPTIX exists to develop products that offer all stakeholders involved in ophthalmic surgery solutions that contribute to a more cost-controlled work environment, as well as improved quality of care.

We are passionate about eliminating cross contamination and reducing risk management whilst increasing efficiency, convenience and reliability. Our innovative solution provides total peace of mind, confidence, and satisfaction.



General information and IP protection

@ HASA OPTIX SRL. All rights reserved. $\ensuremath{\mathbb{B}}/^{\ensuremath{\text{TM}}}$ are trademarks of HASA OPTIX SRL

All original illustrations and text in this publication are fully protected by copyright, no part of it may be reproduced in any form without written permission by the publisher. Size specifications are stated as accurately as possible. Tip drawings are enlarged to show tip detail. Surgical instruments shown in this catalogue are made of stainless steel, unless otherwise indicated.

Why choose HASA OPTIX surgical instruments?

HASA OPTIX is a Belgian medical device company that manufactures a wide range of premium recyclable single-use instruments. Through partnership with HASA OPTIX, Rayner is pleased to offer surgeons a comprehensive portfolio of pre-sterilised instruments and sets which meet high quality material standards.



Premium quality

Premium quality disposable recyclable instruments offering the same constant quality as reusable instruments but with increased convenience and efficiency.



Cost control

HASA OPTIX instruments provide an excellent procurement solution by aiding the precise calculation of surgery costs.



Environmentally friendly

Supplied in recyclable pre-packed pouches relevant to each ophthalmic procedure, helping to reduce waste and environmental impact.



Simplicity

Pre-packed pouches help reduce the number of SKUs, contributing towards a simpler, more efficient repository solution.



Peace of mind

Innovative solutions that provide total peace of mind ensuring the correct surgical instruments are easily accessible.




Infection control

Single use and pre-sterilised instruments ensure all surgeons can be confident with their patient safety of care.

A sustainable and innovative solution with precision performance, that contributes to a better cost-controlled environment and increased quality of patient care.

Instrument Specifications

Capsulorhexis Forcep



23G/25 Gauge - The squeeze handle of this instrument provides delicate, tactile feedback and enables excellent control.


Instrument Specifications

Anatomy Forcep























HANDLES



Flat - This instrument features a flat handle for good manual grip and control of the tip, but may require significant hand, wrist, or forearm movement.



Round - This instrument has ability to rotate the handle that allows for micro control of the tip with minimum movement by rolling between the fingers.



Cross Action - This forcep opens and closes tips by asserting and releasing pressure, offering micro control.



Castroviejo - This forcep is designed for good manual grip and control of the tip, but may require significant hand, wrist, or forearm movement.


Instrument Specifications

Manipulator



Anatomy Speculums

Non-adjustable



MAXIMUM TO MINIMUM WIDTH

Angled

Adjustable

Double-ended

These instruments allow you

to combine two tips of your choice to create a double-ended manipulator/chopper, reducing the number of SKUs needed.



OVERALL LENGTH

BLADE



Open - Allows good access for instruments e.g Phaco handpiece that might otherwise rest on the top bar of a closed blade making surgery difficult



Closed - Good hold on eyelids for procedures where the top bar does not cause a problem e.g trabeculectomy or temporal incision Phaco.



V-shaped - Open wire blades are designed to retain the surgical drape under the eyelids.



Solid - Holds eyelids securely.

MODEL

Size - Different sizes are available (neonatal / infants/adults).

Model - Available in adjustable and non adjustable models.

Approach - Specific models are available for temporal or nasal approaches.

HASA OPTIX's most popular instruments are available in boxes of 10.

Forceps

HA-551001



Toothed Forceps, Straight Shafts With 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm, 1×2 Teeth

HA-551175



Tying Forceps, Curved Shafts With 6mm Tying Platform, Flat Handle, 85mm

HA-551187



Colibri Toothed Forceps, Curved Shafts, Serrated Handle, 0.12mm, 1×2 Teeth

HA-551222



Jeweler Style Forceps, Angled 45°, 5#, Extra Delicate, 0.05×0.01mm Tip

HA-551235



Moorfields Forceps, Blunt Tipped Forceps For Manipulating The Conjunctiva

HA-651000



Capsulorhexis Forceps, Curved Thin Shafts, Flat Cross-Action Handle

HA-651002



Capsulorhexis Forceps, Curved Thin Shafts, Cross-Action, Round Knurled Handle

HA-651071



ICL Cartridge Loading Forceps, Cross-Action

HASA OPTIX's most popular instruments are available in boxes of 10.

Forceps

HA-651072



ICL Forceps, 21g

HA-651143



Capsulorhexis Forceps, Curved Thin Shafts, Cross-Action Ripped Handle

HA-651151



Capsulorhexis Forceps, 45° Angled Thin Shafts, Flat Handle, 85mm, Incsion 1.8mm To 3.2mm

HA-651161



Lenticule Removal Forceps, Serrated Jaws, Angled, Castroviejo Handle

HA-651141



Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle

HA-651147



Capsulorhexis Forceps, 45° Angled Thin Shafts, Round Knurled Handle

HA-651154



Capsulorhexis Forceps, Curved Thin Shafts, Flat Handle, 105mm, Incsion 1.8mm To 3.2mm


HASA OPTIX's most popular instruments are available in boxes of 10.


Micro Vitreoretinal Forceps




Choppers, Manipulators, & Spatulas










ICL Manipulator, Curved, 1.0mm diameter


HASA OPTIX's most popular instruments are available in boxes of 10.

Choppers, Manipulators, & Spatulas




















HASA OPTIX's most popular instruments are available in boxes of 10.

Scissors

HA-553014



Vannas Scissors, Curved, Sharp Pointed Tips, Tip To Pivot Length 10mm

HA-653041



Iris Scissors Hasa Optix, Straight, Sharp

Speculums

HA-555001



Kratz Barraquer Wire Speculum, 13.5mm Open Blades

HA-555007



Wire Speculum, 14mm Open Blades

HA-555023



Lieberman Adjustable Speculum, V-Shaped Wire

HA-555028



Adjustable Speculum, Speculum, 14mm Round Solid Blades

HASA OPTIX's most popular instruments are available in boxes of 10.

Speculums

HA-655035



Speculum Lieberman Fine, Open Blades

HA-655038



Speculum Adjustable Fine Full Blades

HA-655041



Speculum Adjustable Full Rounded Blades

Gauges

HA-558002



Mendez Degree Gauge, Calibrated Every 10 Degrees Form 0°To 180°, External Diameter 19.0mm, lø 12.0mm


HA-SET001		SKU: HA-SET001
051141 €€	651141	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs
OOO Have Con (2)	551001	Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs
	554012	Nucleus Manipulator, Paddle Shaped Tip - 1 pcs
8	653041	Iris Scissors HASA OPTIX, Straight, Sharp - 1 pcs
	655035	Liberman Speculum, Fine, Open Blades - 1 pcs
HA-SET002		SKU: HA-SET002
851141 CE	651141	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs
OOO Mean ON (2)	551001	Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs
	654120	Manipulator Sinskey 0, 7mm Altitude, Tip 0, 2mm - 1 pcs
8	653041	Iris Scissors HASA OPTIX, Straight, Sharp - 1 pcs
	555033	Adjustable Speculum, Lieberman Speculum, 13mm Open Blades - 1 pcs
HA-SET003		SKU: HA-SET003
651141 CE	651141	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs
OOO Mass Only (2)	551001	Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs
	654090	Double Ended Manipulator, Agarwal Phaco Chopper/ Fukasaku Small Pupil Snapper - 1 pcs
8	653041	Iris Scissors HASA OPTIX, Straight, Sharp - 1 pcs
	655035	Liberman Speculum, Fine, Open Blades - 1 pcs

HA-SET004		SKU: HA-SET004
MAA OTIV (2)	651150	Capsulorhexis Forceps, 45° Angled Thin Shafts, Flat handle, 105mm, Incision 1.8mm to 3.2mm - 1 pcs
999 Hanna Carle (2)	551001	Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs
	554012	Nucleus Manipulator, Paddle Shaped Tip - 1 pcs
	555001	Kratz Barraquer Wire Speculum, 13.5mm Open Blades - 1 pcs
HA-SET005		SKU : HA-SET005
051141 €€	651141	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs
Mean Cycle (2)	551001	Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs
OOO Has Que 2	551175	Tying Forceps, Curved Shafts with 6mm Tying Platforms, Flat Handle, 85mm - 1 pcs
	654111	Double-Ended Manipulator Mushroom/Chopper - 1 pcs
	655040	Speculum Lieberman Open Rounded Blades - 1 pcs
HA-SET006		SKU : HA-SET006
ans 2 C	651143	Capsulorhexis Forceps, Curved Thin Shafts, Cross-Action Serrated Handle - 1 pcs
OOO Mana Order (2)	551001	Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs
Hasa Cont. (2)	551122	Notched Forceps, Straight Shafts with Tying Platforms, Flat Handle, 85mm, 0.25mm Notched - 1 pcs
	654111	Double-Ended Manipulator Mushroom/Chopper - 1 pcs
	655035	Liberman Speculum, Fine, Open Blades - 1 pcs

HA-SET007		SKU: HA-SET007
MAGA OFFIX (2)	651150	Capsulorhexis Forceps, 45° Angled Thin Shafts, Flat Handle, 105mm, Incision 1.8mm to 3.2mm - 1 pcs
OOO Mean Colle 2)	551175	Tying Forceps, Curved Shafts with 6mm Tying Platforms, Flat Handle, 85mm - 1 pcs
Oliver Switch	654111	Double-Ended Manipulator Mushroom/Chopper - 1 pcs
	655040	Speculum Lieberman Open Rounded Blades - 1 pcs
HA-SET008		SKU: HA-SET008
051153 C€	651154	Capsulorhexis Forceps, Curved Thin Shafts, Flat Handle, 105mm, Incision 1.8mm to 3.2mm - 1 pcs
OOO Mana Opix (2)	551001	Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs
	554063	Phaco Chopper, Angled Shaft, 1.6mm Tip, Right Hand - 1 pcs
	655040	Speculum Lieberman Open Rounded blades - 1 pcs
HA-SET009		SKU: HA-SET009
Hosp GETS (2)	551122	Notched Forceps, Straight Shafts with Tying Platforms, Flat Handle, 85mm, 0.25mm Notched - 1 pcs
	554063	Phaco Chopper, Angled Shaft, 1.6mm Tip, Right Hand - 1 pcs
	555023	Lieberman Adjustable Speculum, V-Shaped - 1 pcs

HA-SET010		SKU: HA-SET010
051141 €€	651141	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs
OOO Maga Crist (2)	551001	Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs
G G G MASS COTTE (2)	551199	Kelman-McPherson Lens IOL Forceps, Angled Shafts, Flat Handle - 1 pcs
	554012	Nucleus Manipulator, Paddle Shaped Tip - 1 pcs
Married St. and St.	654086	Double-Ended Manipulator Drysdale/Y-Pusher - 1 pcs
	555007	Wire Speculum, 14mm Open Blades - 1 pcs
LIA CETO11		0//1.114.057011
HA-SET011		SKU: HA-SET011
HA-SET011	651141	SKU: HA-SET011 Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs
HA-SET011	651141 551001	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled
HA-SET011		Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs Toothed Forceps, Straight Shafts with 5mm Tying Platforms,
HA-SET011	551001	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs
HA-SET011	551001 554007	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs Lens Manipulator, Y Nucleus Rotator, Angled, Vertical - 1 pcs

HA-SET012		SKU: HA-SET012
65∏47 C€	651147	Capsulorhexis Forceps, 45° Angled Thin Shafts, Round Knurled Handle - 1 pcs
Head Cold (2)	551049	Colibri Toothed Forceps, Curved Shafts With 5mm Tying Platforms, Flat Handle, 0.10mm, 1x2 Teeth - 1 pcs
MASS COTTO (B)	551199	Kelman-McPherson Lens IOL Forceps, Angled Shafts, Flat Handle - 1 pcs
	554012	Nucleus Manipulator, Paddle Shaped Tip - 1 pcs
	654087	Double-Ended Manipulator Cylinder/Y-Pusher - 1 pcs
	555007	Wire Speculum, 14mm Open Blades - 1 pcs
HA-SET013		SKU: HA-SET013
HA-SET013	651147	SKU: HA-SET013 Capsulorhexis Forceps, 45° Angled Thin Shafts, Round Knurled Handle - 1 pcs
HA-SET013	651147 551049	Capsulorhexis Forceps, 45° Angled Thin Shafts,
HA-SET013		Capsulorhexis Forceps, 45° Angled Thin Shafts, Round Knurled Handle - 1 pcs Colibri Toothed Forceps, Curved Shafts With 5mm Tying
HA-SET013	551049	Capsulorhexis Forceps, 45° Angled Thin Shafts, Round Knurled Handle - 1 pcs Colibri Toothed Forceps, Curved Shafts With 5mm Tying Platforms, Flat Handle, 0.10mm, 1x2 Teeth - 1 pcs
HA-SET013	551049 554007	Capsulorhexis Forceps, 45° Angled Thin Shafts, Round Knurled Handle - 1 pcs Colibri Toothed Forceps, Curved Shafts With 5mm Tying Platforms, Flat Handle, 0.10mm, 1x2 Teeth - 1 pcs Lens Manipulator, Y Nucleus Rotator, Angled, Vertical - 1 pcs



HA-SET016		SKU: HA-SET016
MAGA OFFIN (g)	651151	Capsulorhexis Forceps, 45° Angled Thin Shafts,Flat Handle, 85mm, Incision 1.8mm to 3.2mm - 1 pcs
Mana Color (g)	551025	Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Castroviejo Handle, 0.12mm, 1x2 Teeth - 1 pcs
I B B B HARA GITTE (1)	551199	Kelman-McPherson Lens IOL Forceps, Angled Shafts, Flat Handle - 1 pcs
	554008	Lens Manipulator V-Shaped Tip, Angled - 1 pcs
	554015	Nucleus Manipulator, Connor Wand, Angled, Round - 1 pcs
	654086	Double-Ended Manipulator Drysdale/Y-Pusher - 1 pcs
	553039	Eye Scissors, Straight, Blunt Ended, 90mm - 1 pcs
	555023	Lieberman Adjustable Speculum, V-Shaped Wire - 1 pcs
HA-SET017		SKU: HA-SET017
HA-SET017		
HA-SET017	651141	SKU: HA-SET017 Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs
HA-SET017	651141 551001	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled
HA-SET017		Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs Toothed Forceps, Straight Shafts with 5mm Tying Platforms,
651/41 C¢	551001	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs Colibri Toothed Forceps, Curved Shafts, Serrated Handle,
651/41 C¢	551001 551187	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs Colibri Toothed Forceps, Curved Shafts, Serrated Handle, 0.12mm 1x2 teeth - 1 pcs Akahoshi Prechopper Forceps, Cross-Action, 10mm Angled
651/41 C¢	551001 551187 551206	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs Colibri Toothed Forceps, Curved Shafts, Serrated Handle, 0.12mm 1x2 teeth - 1 pcs Akahoshi Prechopper Forceps, Cross-Action, 10mm Angled Shafts - 1 pcs

HA-SET019		SKU: HA-SET019
05141 C€	651141	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs
OOO Maga Cork (2)	551001	Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs
Aba a Critic (2)	551175	Tying Forceps, Curved Shafts with 6mm Tying Platforms, Flat Handle, 85mm - 1 pcs
	653042	Iris Scissors HASA OPTIX, Curved, Sharp - 1 pcs
	655035	Liberman Speculum, Fine, Open Blades - 1 pcs
HA-SET022		SKU: HA-SET022
MAGA OTTA (g)	651150	Capsulorhexis Forceps, 45° Angled Thin Shafts, Flat Handle, 105mm, Incision 1.8mm to 3.2mm - 1 pcs
MAG DITE (2)	551187	Colibri Toothed Forceps, Curved Shafts, Serrated Handle, 0.12mm 1x2 teeth - 1 pcs
	654117	Double-Ended Manipulator Drysdale/Mushroom - 1 pcs
	654120	Manipulator Sinskey 0, 7mm Altitude, Tip 0, 2mm - 1 pcs
8	653041	Iris Scissors HASA OPTIX, Straight, Sharp - 1 pcs
	555017	Wire Speculum, 14mm Open Blades - 1 pcs
0	558002	Mendez Degree Gauge, Calibrated every 10 degrees from 0° to 180° - 1 pcs
HA-SET024		SKU: HA-SET024
		ONO . HA 3E1024
	651075	IOL Cutting Forceps, 23G
	653059	IOL Cutting Scissors 19G



HA-SET033		SKU: HA-SET033
JACA 1995 P	651161	Lenticule Removal Forceps, Serrated Jaws, Angled
	654108	Lenticule Separator/Elevator
O O O Mass Code (2)	551175	Tying Forceps, Curved Shafts with 6mm Tying Platforms, Flat Handle, 85mm - 1 pcs
	555031	Lindstrom-Chu Aspirating Speculum
LIA OFTOOA		
HA-SET034		SKU: HA-SET034
Maga Ciris (2)	551001	Toothed Forceps, Straight Shafts with 5mm Tying Platforms, Flat Handle, 85mm, 0.12mm 1x2 Teeth - 1 pcs
051141 C€	651141	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1 pcs
	555007	Wire Speculum, 14mm Open Blades - 1 pcs
	654112	Manipulator, Cylinder
HA-SET035		SKU: HA-SET035
OOO Hage COV. (2)	551001	Toothed Forceps Straight shafts with 5mm tying platforms flat handle 0.12mm teeth
051141 CC	651141	Capsulorhexis Forceps, Curved Thin Shafts, Round Knurled Handle - 1pcs
	555007	Wire Speculum, 14mm Open Blades - 1 pcs
	654096	Double-Ended Manipulator, Drysdale/Koch Chopper



Join us to simplify ophthalmic surgery



Please visit **rayner.com/hasa** or contact your local Rayner representative for more information

The HASA Optix App is available on Google Play & App Store

CLICK BELOW TO DOWNLOAD



